


Product Datasheet

DNAJC13 Antibody (orb1239876)

Description

DNAJC13 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	DNAJC13 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human DNAJC13. The immunogen is located within amino acids 1940 - 1990 of DNAJC13.
Target	DNAJC13
Preservatives	DNAJC13 antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	DNAJC13 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	DNAJC13 antibody can be used for detection of DNAJC13 by Western blot at 1 µg/mL. Antibody can also be used for immunohistochemistry starting at 5 µg/mL. Antibody validated: Western Blot in mouse samples and Immunohistochemistry in human samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 249 kDa Observed: 249 kDa
Uniprot ID	O75165
NCBI	NP_056083
Dilution Range	DNAJC13 antibody can be used for detection of DNAJC13 by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunohistochemistry starting at 5 µg/mL. Antibody validated: Western Blot in human samples and Immunohistochemistry in human samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



WB Validation of
DNAJC13 with Mouse
Samp...



Immunohistochemistry
of DNAJC13 in
human...