



Product Datasheet

DDDDK-tag Antibody [1D1B12] (orb1239896)

Description	DDDDK-tag Antibody [1D1B12]
Species/Host	Mouse
Reactivity	Tag
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	DDDDK-tag antibody was raised against a synthetic peptide containing DYKDDDDK.
Preservatives	DDDDK-tag antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	DDDDK-tag antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	DDDDK-tag antibody can be used for detection of DDDDK-tag epitope-containing recombinant proteins by Western blot at 0.5 - 1 µg/ml.
Isotype	IgG
Clonality	Monoclonal
Dilution Range	DDDDK-tag antibody can be used for detection of DDDDK-tag epitope-containing recombinant proteins by Western blot at 0.5 - 1 µg/ml.
Expiration Date	12 months from date of receipt.



Western blot analysis of (1) 1000ng, (2)...



Western blot analysis of (1) 20 ng of DD...