

Product Datasheet

PPARGC1A Antibody (orb1239904)

Description

PPARGC1A Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

PPARGC1A antibody was raised against a 17 amino acid peptide near the amino terminus of human PPARGC1A. The immunogen is located within amino acids 200 - 250 of PPARGC1A.

Target

PPARGC1A

Preservatives

PPARGC1A antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

PPARGC1A antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

WB: 2 µg/mL; IHC: 1-2 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunohistochemistry in human, mouse and rat samples; Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 91kD Observed: 100 kD

Uniprot ID

[Q9UBK2](#)

NCBI


[NP_037393](#)

Dilution Range


WB: 2 µg/mL; IHC: 1-2 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunohistochemistry in human, mouse and rat samples; Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



WB
Validation in
Human Cell
Lines. Loadi...



WB
Validation in
Human Lung.
Loading: 15...



WB
Validation in
Mouse Skin.
Loading: 15...