

Product Datasheet

CNRIP1 Antibody (orb1239963)

Description

CNRIP1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

CNRIP1 antibody was raised against a 15 amino acid peptide near the center of human CNRIP1. The immunogen is located within amino acids 50 - 100 of CNRIP1.

Target

CNRIP1

Preservatives

CNRIP1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

CNRIP1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

CNRIP1 antibody can be used for detection of CNRIP1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 18 kDa Observed: 16 kDa

Uniprot ID

[Q96F85](#)

NCBI


[NP_056278](#)

Dilution Range


CNRIP1 antibody can be used for detection of CNRIP1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western blot analysis of CNRIP1 in human...



Immunohistochemistry of CNRIP1 in rat br...



Immunofluorescence of CNRIP1 in rat brai...