


Product Datasheet


SARS-CoV-2 (COVID-19) Nucleocapsid Antibody (orb1239972)

Description


SARS-CoV-2 (COVID-19) Nucleocapsid



Immunohistochemistry
Validation of SARS-...



B) or SARS-CoV-2
Spike S1 antibodies
(or...



E) in adjacent sections
of autopsy lung ...

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Anti-SARS-CoV-2 (COVID-19) Nucleocapsid antibody (orb1239972) was raised against a peptide corresponding to 17 amino acids near the center of SARS-CoV-2 (COVID-19) Nucleocapsid protein. The immunogen is located within 230-280 amino acids of SARS-CoV-2 (COVID-19) Nucleocapsid protein.

Target

N

Preservatives

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SARS-CoV-2 (COVID-19) Nucleocapsid antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

NCBI

[QHD43423](#)

Dilution Range

WB: 0.5-2 µg/mL SARS-CoV-2 (COVID-19) Nucleocapsid antibody can be used for the detection of SARS-CoV-2 (COVID-19) Nucleocapsid protein in ELISA. It will detect 4 ng of free peptide at 1 µg/mL. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.