

Product Datasheet

CITED2 Antibody (orb1239973)

Description
CITED2 Antibody
Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, WB

Applications
Immunogen

CITED2 antibody was raised against an 18 amino acid synthetic peptide near the amino terminus of human CITED2. The immunogen is located within the first 50 amino acids of CITED2.

Target

CITED2

Preservatives

CITED2 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

CITED2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

CITED2 antibody can be used for detection of CITED2 by Western blot at 1 - 2 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality


Polyclonal

MW


30 kDa

Uniprot ID
Q99967
NCBI
NP_001161861
Dilution Range

CITED2 antibody can be used for detection of CITED2 by Western blot at 1 - 2 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.



Western blot analysis of CITED2 in Jurka...



Immunofluorescence of CITED2 in Jurkat C...