


Product Datasheet

CEACAM5 Antibody (orb1239991)


Description

CEACAM5 Antibody

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	CEACAM5 antibody was raised against a 15 amino acid peptide near the center of human CEACAM5. The immunogen is located within amino acids 310 - 360 of CEACAM5 .
Target	CEACAM5
Preservatives	CEACAM5 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	CEACAM5 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	CEACAM5 antibody can be used for detection of CEACAM5 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in rat samples and Immunofluorescence in human samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 77 kDa Observed: 80kDa
Uniprot ID	P06731
NCBI	NP_004354
Dilution Range	CEACAM5 antibody can be used for detection of CEACAM5 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in rat samples and Immunofluorescence in human samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Western blot analysis of CEACAM5 in rat ...



Immunofluorescence of CEACAM5 in human l...