




Product Datasheet

Diablo Antibody (orb1240005)

Description	Diablo Antibody	
Species/Host	Rabbit	 Western Blot Validation in (A and B) Mou...
Conjugation	Unconjugated	
Tested Applications	ELISA, IF, IHC-P, IP, WB	
Immunogen	Anti-Smac antibody (orb1240005) was raised against a peptide corresponding to 16 amino acids near the carboxy terminus of murine Smac. The immunogen is located within the last 50 amino acids of Smac.	
Target	Diablo	 Independent Antibody Validation (IAV) vi...
Preservatives	Smac Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/ml	
Storage	Smac antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	 Western Blot Validation in Human, Mouse ...
Note	For research use only	
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 27kD Observed: 20 kD	
Uniprot ID	Q9JIQ3	
NCBI	AF203914	
Dilution Range	WB: 1 µg/mL; IF: 10 µg/mL; IHC: 2 µg/mL. Antibody validated: Western Blot in mouse and rat samples; Immunofluorescence and Immunohistochemistry in mouse samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	