

Product Datasheet

CCNT1 Antibody (orb1240019)

Description

CCNT1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

Rabbit polyclonal CCNT1 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human CCNT1. The immunogen is located within the last 50 amino acids of CCNT1.

Target

CCNT1

Preservatives

CCNT1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

CCNT1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

CCNT1 antibody can be used for detection of CCNT1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 80 kDa

Uniprot ID

[O60563](#)

NCBI


[NP_001231](#)

Dilution Range


CCNT1 antibody can be used for detection of CCNT1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western blot analysis of CCNT1 in rat br...



Immunohistochemistry of CCNT1 in rat bra...



Immunofluorescence of CCNT1 in rat brain...