

Product Datasheet

CASP10 Antibody (orb1240057)

Description

CASP10 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, WB

Applications

Immunogen

Caspase-10 antibody was raised against a 17 amino acid peptide near the carboxy terminus of human Caspase-10. The immunogen is located within amino acids 240 - 290 of Caspase-10.

Target

CASP10

Preservatives

Caspase-10 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/ml

Storage

Caspase-10 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

Caspase-10 antibody can be used for detection of caspase-10 by Western blot 1 µg/mL. A 59 kDa band can be detected. Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

[Q92851](#)


NCBI

[AAD28402](#)


Dilution Range

Caspase-10 antibody can be used for detection of caspase-10 by Western blot 1 µg/mL. A 59 kDa band can be detected. Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.


H J A K



Western blot analysis of caspase-10 in H...



Immunofluorescence of Casp10 in Hela cel...



Immunocytochemistry of Casp10 in Hela ce...