

Product Datasheet

BTLA Antibody (orb1240080)

Description

BTLA Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

BTLA antibody was raised against a 17 amino acid peptide near the center of human BTLA. The immunogen is located within amino acids 170 - 220 of BTLA .

Target

BTLA

Preservatives

BTLA Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

BTLA antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

BTLA antibody can be used for detection of BTLA by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 32 kDa Observed: 40kDa

Uniprot ID

[Q7Z6A9](#)

NCBI

[NP_861445](#)


Dilution Range

BTLA antibody can be used for detection of BTLA by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.


Expiration Date

12 months from date of receipt.


A B



Western blot analysis of BTLA in Jurkat ...



Immunohistochemistry of BTLA in human sm...



Immunofluorescence of BTLA in human smal...