

## Product Datasheet

### BAK1 Antibody (orb1240143)

## Description

BAK1 Antibody

### Species/Host

Rabbit

### Conjugation

Unconjugated

### Tested Applications

ELISA, ICC, WB

### Immunogen

Bak antibody was raised against a peptide corresponding to 13 amino acids near the amino-terminus of human Bak. The immunogen is located within amino acids 30 - 80 of Bak.

### Target

BAK1

### Preservatives

Bak Antibody is supplied in PBS containing 0.02% sodium azide.

### Form/Appearance

Liquid

### Concentration

1 mg/mL

### Storage

Bak antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

Bak antibody can be used for detection of Bak by Western blot at 1 to 2 µg/mL. Antibody can also be used for immunocytochemistry starting at 2 µg/mL. Antibody validated: Western Blot in mouse samples and Immunocytochemistry in mouse samples. All other applications and species not yet tested.

### Isotype

IgG

### Clonality

Polyclonal

### Uniprot ID

[Q16611](#)

### NCBI


[Q16611](#)

### Dilution Range


Bak antibody can be used for detection of Bak by Western blot at 1 to 2 µg/mL. Antibody can also be used for immunocytochemistry starting at 2 µg/mL. Antibody validated: Western Blot in mouse samples and Immunocytochemistry in mouse samples. All other applications and species not yet tested.

### Expiration Date

12 months from date of receipt.



Western blot analysis of Bak in L1210 ce...



Immunocytochemistry of Bak in L1210 cell...