

Product Datasheet

HAVCR1 Antibody (orb1240149)

Description

HAVCR1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

Anti-TIM-1 antibody (orb1240149) was raised against a peptide corresponding to 16 amino acids near the amino terminus of human TIM-1. The immunogen is located within amino acids 50 - 100 of TIM-1.

Target

HAVCR1

Preservatives

TIM-1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

TIM-1 antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 44kD Observed: 60kD (Post-modification: 4 N-linked glycosylation)

Uniprot ID

[Q96D42](#)

NCBI


[NP_036338](#)

Dilution Range


WB: 1 - 8 µg/mL (overnight incubation at 4° C); IHC-P: 10 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human and mouse samples; Immunohistochemistry and Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western Blot
Validation in
Human Cell
Li...



Independent
Antibody
Validation
(IAV) vi...



Validation
with TIM-1
siRNA
Knockdown
in...