

Product Datasheet

APOBEC3B Antibody (orb1240220)

Description

APOBEC3B Antibody

KDa
95 -

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

APOBEC3B antibody was raised against a 14 amino acid peptide near the amino terminus of human APOBEC3B. The immunogen is located within amino acids 20 - 70 of APOBEC3B.

Target

APOBEC3B

Preservatives

APOBEC3B Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

APOBEC3B antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

APOBEC3B antibody can be used for detection of APOBEC3B by Western blot at 1 - 2 µg/ml. Antibody can also be used for Immunohistochemistry at 5 µg/mL. For Immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 42 kDa Observed: 48 kDa

Uniprot ID



[Q9UH17](#)

NCBI


[NP_004891](#)

Dilution Range

APOBEC3B antibody can be used for detection of APOBEC3B by Western blot at 1 - 2 µg/ml. Antibody can also be used for Immunohistochemistry at 5 µg/mL. For Immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All



Immunohistochemistry of APOBEC3B in human...



Immunofluorescence of APOBEC3B in human ...