

Product Datasheet

ALSFTD Antibody (orb1240252)

Description

ALSFTD Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, WB

Applications

Immunogen

ALSFTD antibody was raised against a 19 amino acid peptide from near the carboxy terminus of human ALSFTD. The immunogen is located within amino acids 400 - 450 of ALSFTD.

Target

ALSFTD

Preservatives

ALSFTD Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

ALSFTD antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

ALSFTD antibody can be used for the detection of ALSFTD by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunocytochemistry at 10 µg/mL. Antibody validated: Western Blot in human samples and Immunocytochemistry in mouse samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 53 kDa Observed: 52 kDa

Uniprot ID

[Q96LT7](#)

NCBI


[NP_060795](#)

Dilution Range


ALSFTD antibody can be used for the detection of ALSFTD by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunocytochemistry at 10 µg/mL. Antibody validated: Western Blot in human samples and Immunocytochemistry in mouse samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of ALSFTD in HeLa ...



Immunocytochemistry of ALSFTD in A-20 ce...