

Product Datasheet

TNFSF12 Antibody (orb1240271)

Description

TNFSF12 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

ELISA, IF, IHC-P, WB

Immunogen

Anti-TWEAK antibody (orb1240271) was raised against recombinant human TWEAK protein.

Target

TNFSF12

Preservatives

TWEAK Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

TWEAK antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 27 and 36 kDa Observed: 27 and 36 kDa

Uniprot ID

O43508

NCBI


NP_003800

Dilution Range


WB: 0.5-4 µg/mL; IHC: 5-10 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human and mouse samples; Immunohistochemistry in human and mouse samples; Immunofluorescence in mouse samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western Blot Validation in HeLa Cell Lys...



Western Blot Validation with Recombinant...



Western Blot Validation in Mouse Pancrea...