




Product Datasheet

UNG Antibody (orb1240273)

Description	UNG Antibody	
Species/Host	Rabbit	 Western Blot Validation in Mouse C2C12 C...
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, IF, IHC-P, WB	
Immunogen	Anti-UNG1 antibody (orb1240273) was raised against a peptide corresponding to 13 amino acids near the carboxy terminus of human UNG1. The immunogen is located within the last 50 amino acids of UNG1.	
Target	UNG	 Independent Antibody Validation (IAV) vi...
Preservatives	UNG1 Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	UNG1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	 Western Blot Validation in HeLa Cell Lys...
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 34kD Observed: 34 kD	
Uniprot ID	P13051	
NCBI	NP_003353	
Dilution Range	WB: 0.5-4 µg/mL; IHC: 2 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human and mouse samples; Immunohistochemistry in human samples; Immunofluorescence in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	