

## Product Datasheet

### AGR2 Antibody (orb1240275)

**Description**
**AGR2 Antibody**
**Species/Host**

Rabbit

**Reactivity**

Human, Mouse

**Conjugation**

Unconjugated

**Tested**

ELISA, IF, IHC-P, WB

**Applications**
**Immunogen**

Rabbit polyclonal AGR2 antibody was raised against a 17 amino acid peptide near the amino terminus of human AGR2. The immunogen is located within the first 50 amino acids of AGR2.

**Target**

AGR2

**Preservatives**

AGR2 Antibody is supplied in PBS containing 0.02% sodium azide.

**Form/Appearance**

Liquid

**Concentration**

1 mg/mL

**Storage**

AGR2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

**Note**

For research use only

**Application notes**

AGR2 antibody can be used for detection of AGR2 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**


19 kDa

**Uniprot ID**
[O95994](https://www.uniprot.org/entry/O95994)
**NCBI**
[NP\\_006399](https://www.ncbi.nlm.nih.gov/nuccore/NP_006399)
**Dilution Range**


AGR2 antibody can be used for detection of AGR2 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

**Expiration Date**


12 months from date of receipt.

**A B**


Western blot analysis of AGR2 in Hela ce...



Immunohistochemistry of AGR2 in human sm...



Immunofluorescence of AGR2 in human smal...