



Product Datasheet

EZH2 antibody (orb12403)


Description	Rabbit polyclonal antibody to EZH2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human EZH2/KMT6(NP_001190176.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	78-110kDa
Uniprot ID	Q15910
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:20-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric cancer tissue.



Immunohistochemical staining of human brain tissue.



Immunohistochemical staining of rat brain tissue.