

Product Datasheet

ACE2 Antibody (HRP) (orb1240300)

Description

ACE2 Antibody (HRP)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

HRP

Tested

ELISA, IHC

Applications

Immunogen

Anti-ACE2 antibody (3217) was raised against a peptide corresponding to 18 amino acids near the carboxy terminus of human ACE2. The immunogen is located within the last 50 amino acids of ACE2.

Target

ACE2

Preservatives

ACE2 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

ACE2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

IHC: 2 µg/mL Antibody validated:
Immunohistochemistry in human samples. All other applications and species not yet tested

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 93kD
Observed: 130 kD (7 N-linked glycosylation)

Uniprot ID

[Q9BYF1](#)

NCBI


[NP_068576](#)

Dilution Range


IHC: 2 µg/mL Antibody validated:
Immunohistochemistry in human samples. All other applications and species not yet tested

Expiration Date

12 months from date of receipt.



ELISA Validation.
Coating Antigen:
immun...



Immunohistochemistry
Validation of ACE-2...