


Product Datasheet

**SARS-CoV-2 (COVID-19) S protein RBD antibody
[F3C10] (orb1240313)**

Description

SARS-CoV-2 (COVID-19) S protein RBD

PROBIO
TEST TABLE



orb1240313
binds with
RBD
Variants

Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Recombinant SARS-CoV-2 S Protein RBD. This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant SARS-CoV-2(2019-nCoV) RBD protein.
Target	S
Preservatives	PBS, pH 7.4, containing 0.05% proclin300, 50% glycerol.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage.
Note	For research use only
Application notes	WB:1:1000 - 1:3000; ELISA:1:1000-1:3000; ICC:1:1000.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P0DTC2
NCBI	P0DTC2
Dilution Range	WB:1:1000 - 1:3000; ELISA:1:1000-1:3000; ICC:1:1000.
Expiration Date	12 months from date of receipt.