

## Product Datasheet

### RHOT1 Antibody (orb1240334)

## Description

## RHOT1 Antibody

### Species/Host

Rabbit

### Reactivity

Canine, Human, Mouse

### Conjugation

Unconjugated

### Tested

ELISA, IHC, WB

### Applications

### Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RHOT1.

### Target

RHOT1

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store RHOT1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

### Note

For research use only

### Application notes

RHOT1 antibody can be used for detection of RHOT1 by ELISA at 1:62500. RHOT1 antibody can be used for detection of RHOT1 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

71 kDa, 76 kDa, 80 kDa

### Uniprot ID

[Q8IXI2](#)

### NCBI


[NP\\_060777](#)

### Dilution Range


RHOT1 antibody can be used for detection of RHOT1 by ELISA at 1:62500. RHOT1 antibody can be used for detection of RHOT1 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Muscle at ...



See protocol 2204 and 2205 for more information

Antibody used in WB on mouse brain and h...



Antibody used in WB on Human Fetal Lung ...