

Product Datasheet

TM9SF1 Antibody (orb1240349)

Description

TM9SF1 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TM9SF1.

Target

TM9SF1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TM9SF1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TM9SF1 antibody can be used for detection of TM9SF1 by ELISA at 1:312500. TM9SF1 antibody can be used for detection of TM9SF1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

67 kDa, 55 kDa

Uniprot ID

[O15321](#)

NCBI


[NP_006396](#)

Dilution Range

TM9SF1 antibody can be used for detection of TM9SF1 by ELISA at 1:312500. TM9SF1 antibody can be used for detection of TM9SF1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in WB on Human Jurkat 1.25...