

Product Datasheet

PDPN Antibody (orb1240356)

Description

PDPN Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PDPN.

Target

PDPN

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PDPN antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PDPN antibody can be used for detection of PDPN by ELISA at 1:62500. PDPN antibody can be used for detection of PDPN by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

25 kDa, 12 kDa, 12 kDa, 25 kDa

Uniprot ID

[Q86YL7](#)

NCBI


[NP_006465](#)

Dilution Range


PDPN antibody can be used for detection of PDPN by ELISA at 1:62500. PDPN antibody can be used for detection of PDPN by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Lung.



Antibody used in WB on Human Jurkat 0.2-...