

Product Datasheet

CHAC2 Antibody (orb1240362)

Description

CHAC2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Placenta at...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CHAC2.

Target

CHAC2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CHAC2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CHAC2 antibody can be used for detection of CHAC2 by ELISA at 1:12500. CHAC2 antibody can be used for detection of CHAC2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

21 kDa

Uniprot ID

[Q8WUX2](#)

NCBI

[NP_001008708](#)

Dilution Range

CHAC2 antibody can be used for detection of CHAC2 by ELISA at 1:12500. CHAC2 antibody can be used for detection of CHAC2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.