

Product Datasheet

CHIA Antibody (orb1240409)

Description

CHIA Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CHIA.

Target

CHIA

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CHIA antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CHIA antibody can be used for detection of CHIA by ELISA at 1:62500. CHIA antibody can be used for detection of CHIA by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

40 kDa

Uniprot ID

[Q9BZP6](#)

NCBI

[NP_068569](#)

Dilution Range

CHIA antibody can be used for detection of CHIA by ELISA at 1:62500. CHIA antibody can be used for detection of CHIA by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

91kDa

41kDa

31kDa

21kDa

Antibody
used in WB
on
Transfected
293T ...

CHIA

91kDa

65kDa

40kDa

31kDa

22kDa

Rabbit Anti-CHIA
Sample Type: Human Fetal Heart
Antibody Concentration: 1µg/mL

Antibody
used in WB
on Hum. Fetal
Heart ...

CHIA

91kDa

65kDa

40kDa

31kDa

22kDa

Rabbit Anti-CHIA
Sample Type: Human Fetal Lung
Antibody Concentration: 1µg/mL

Antibody
used in WB
on Hum. Fetal
Lung a...