

Product Datasheet


PLAG1 Antibody (orb1240415)

Description

PLAG1 Antibody

Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PLAG1.
Target	PLAG1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PLAG1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	PLAG1 antibody can be used for detection of PLAG1 by ELISA at 1:312500. PLAG1 antibody can be used for detection of PLAG1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	Q6DJT9
NCBI	NP_002646
Dilution Range	PLAG1 antibody can be used for detection of PLAG1 by ELISA at 1:312500. PLAG1 antibody can be used for detection of PLAG1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.


PLAG1



Antibody used in WB on Hum. Fetal Lung a...



Antibody used in WB on Human OVCAR-3 at ...



Antibody used in WB on Hum. Fetal Muscle...