

Product Datasheet

TAF7 Antibody (orb1240448)

Description

TAF7 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TAF7.

Target

TAF7

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TAF7 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TAF7 antibody can be used for detection of TAF7 by ELISA at 1:1562500. TAF7 antibody can be used for detection of TAF7 by western blot at 5.0-8.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

40 kDa

Uniprot ID

[Q15545](#)

NCBI


[NP_005633](#)

Dilution Range


TAF7 antibody can be used for detection of TAF7 by ELISA at 1:1562500. TAF7 antibody can be used for detection of TAF7 by western blot at 5.0-8.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in IHC on Human Lung.



Antibody used in WB on Human Jurkat 5.0-...