

Product Datasheet

SNAPC2 Antibody (orb1240451)

Description

SNAPC2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SNAPC2.

Target

SNAPC2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SNAPC2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SNAPC2 antibody can be used for detection of SNAPC2 by ELISA at 1:1562500. SNAPC2 antibody can be used for detection of SNAPC2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

36 kDa

Uniprot ID

[Q13487](#)

NCBI


[NP_003074](#)

Dilution Range

SNAPC2 antibody can be used for detection of SNAPC2 by ELISA at 1:1562500. SNAPC2 antibody can be used for detection of SNAPC2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in WB on Human Jurkat 1.25...