

Product Datasheet

POU3F4 Antibody (orb1240469)

Description

POU3F4 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
HepG2 at
0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human POU3F4.

Target

POU3F4

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store POU3F4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

POU3F4 antibody can be used for detection of POU3F4 by ELISA at 1:1562500. POU3F4 antibody can be used for detection of POU3F4 by western blot at 0.3-0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

39 kDa

Uniprot ID

[P49335](#)

NCBI

[NP_000298](#)

Dilution Range

POU3F4 antibody can be used for detection of POU3F4 by ELISA at 1:1562500. POU3F4 antibody can be used for detection of POU3F4 by western blot at 0.3-0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.