


Product Datasheet

IRF8 Antibody (orb1240488)


Description

IRF8 Antibody


Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human IRF8.
Target	IRF8
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store IRF8 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	IRF8 antibody can be used for detection of IRF8 by ELISA at 1:62500. IRF8 antibody can be used for detection of IRF8 by western blot at 0.5-2.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	48 kDa
Uniprot ID	Q02556
NCBI	NP_002154
Dilution Range	IRF8 antibody can be used for detection of IRF8 by ELISA at 1:62500. IRF8 antibody can be used for detection of IRF8 by western blot at 0.5-2.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Intestine.



Antibody used in IHC on Human Liver.



Antibody used in WB on Human Jurkat 0.2-...