

Product Datasheet

CREB5 Antibody (orb1240516)

Description

CREB5 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
HepG2 at
1....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CREB5.

Target

CREB5

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CREB5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CREB5 antibody can be used for detection of CREB5 by ELISA at 1:62500. CREB5 antibody can be used for detection of CREB5 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

53 kDa, 56 kDa, 57 kDa

Uniprot ID

[Q02930](#)

NCBI

[NP_878902](#)

Dilution Range

CREB5 antibody can be used for detection of CREB5 by ELISA at 1:62500. CREB5 antibody can be used for detection of CREB5 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.