


Product Datasheet

SNAP25 Antibody (orb1240546)


Description

SNAP25 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse SNAP25.
Target	SNAP25
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SNAP25 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SNAP25 antibody can be used for detection of SNAP25 by ELISA at 1:312500. SNAP25 antibody can be used for detection of SNAP25 by western blot at 1.0μg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	23 kDa
Uniprot ID	P60879
NCBI	NP_035558
Dilution Range	SNAP25 antibody can be used for detection of SNAP25 by ELISA at 1:312500. SNAP25 antibody can be used for detection of SNAP25 by western blot at 1.0μg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Frontal Cor...



Antibody used in WB on Human NT-2 and Mo...