

Product Datasheet

HTR1A Antibody (orb1240566)

Description

HTR1A Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human HTR1A.

Target

HTR1A

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store HTR1A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

HTR1A antibody can be used for detection of HTR1A by ELISA at 1:62500. HTR1A antibody can be used for detection of HTR1A by western blot at 0.0625 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

46 kDa

Uniprot ID

[P08908](#)

NCBI


[NP_000515](#)

Dilution Range


HTR1A antibody can be used for detection of HTR1A by ELISA at 1:62500. HTR1A antibody can be used for detection of HTR1A by western blot at 0.0625 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human myometrial.



Antibody used in WB on Human Jurkat 1 ug...



Antibody used in WB on Jurkat 1.0 ug/ml ...