

Product Datasheet

CHRNE Antibody (orb1240577)

Description		7500
Species/Host	Rabbit	4500
Reactivity	Human	3500
Conjugation	Unconjugated	2500
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CHRNE.	Antibody used in WB on Human HepG2 at 1....
Target	CHRNE	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CHRNE antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	CHRNE antibody can be used for detection of CHRNE by ELISA at 1:375000. CHRNE antibody can be used for detection of CHRNE by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	55 kDa	
Uniprot ID	Q04844	
NCBI	NP_000071	
Dilution Range	CHRNE antibody can be used for detection of CHRNE by ELISA at 1:375000. CHRNE antibody can be used for detection of CHRNE by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	