


Product Datasheet

CHRND Antibody (orb1240579)

Description

CHRND Antibody

Species/Host	Rabbit
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CHRND.
Target	CHRND
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CHRND antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	CHRND antibody can be used for detection of CHRND by ELISA at 1:62500. CHRND antibody can be used for detection of CHRND by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	59 kDa
Uniprot ID	Q07001
NCBI	NP_000742
Dilution Range	CHRND antibody can be used for detection of CHRND by ELISA at 1:62500. CHRND antibody can be used for detection of CHRND by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in WB on Human HepG2 at 0....