

Product Datasheet

CHRNA9 Antibody (orb1240583)

Description

CHRNA9 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CHRNA9.

Target

CHRNA9

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CHRNA9 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CHRNA9 antibody can be used for detection of CHRNA9 by ELISA at 1:312500. CHRNA9 antibody can be used for detection of CHRNA9 by western blot at 0.0625 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

55 kDa

Uniprot ID

[Q9UGM1](#)

NCBI

[NP_060051](#)

Dilution Range

CHRNA9 antibody can be used for detection of CHRNA9 by ELISA at 1:312500. CHRNA9 antibody can be used for detection of CHRNA9 by western blot at 0.0625 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human Jurkat 0.06...