

Product Datasheet

CCNH Antibody (orb1240636)

Description

CCNH Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CCNH.

Target

CCNH

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CCNH antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CCNH antibody can be used for detection of CCNH by ELISA at 1:75000. CCNH antibody can be used for detection of CCNH by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

38 kDa

Uniprot ID

[P51946](#)

NCBI


[NP_001230](#)

Dilution Range


CCNH antibody can be used for detection of CCNH by ELISA at 1:75000. CCNH antibody can be used for detection of CCNH by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Lung.



Antibody used in WB on Human Jurkat 1 ug...