

Product Datasheet

TNFRSF18 Antibody (orb1240646)

Description	TNFRSF18 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human TNFRSF18.	Antibody used in WB on Human HeLa at 1 u...
Target	TNFRSF18	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store TNFRSF18 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	TNFRSF18 antibody can be used for detection of TNFRSF18 by ELISA at 1:312500. TNFRSF18 antibody can be used for detection of TNFRSF18 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	23 kDa	
Uniprot ID	B1AME3	
NCBI	NP_683700	
Dilution Range	TNFRSF18 antibody can be used for detection of TNFRSF18 by ELISA at 1:312500. TNFRSF18 antibody can be used for detection of TNFRSF18 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	