




Product Datasheet

TSKS Antibody (orb1240701)

Description	TSKS Antibody	
Species/Host	Rabbit	 
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TSKS.	Antibody used in IHC on Human Testis at ...
Target	TSKS	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store TSKS antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Human HepG2 at 1 ...
Note	For research use only	
Application notes	TSKS antibody can be used for detection of TSKS by ELISA at 1:1562500. TSKS antibody can be used for detection of TSKS by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	65 kDa	
Uniprot ID	Q9UJT2	Antibody used in IHC on Human Testis.
NCBI	NP_068379	
Dilution Range	TSKS antibody can be used for detection of TSKS by ELISA at 1:1562500. TSKS antibody can be used for detection of TSKS by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	