

Product Datasheet

SAMD4A Antibody (orb1240709)

Description

SAMD4A Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SAMD4A.

Target

SAMD4A

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SAMD4A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SAMD4A antibody can be used for detection of SAMD4A by ELISA at 1:1562500. SAMD4A antibody can be used for detection of SAMD4A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

79 kDa

Uniprot ID

[Q9UPU9](#)

NCBI

[NP_056404](#)

Dilution Range

SAMD4A antibody can be used for detection of SAMD4A by ELISA at 1:1562500. SAMD4A antibody can be used for detection of SAMD4A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90 kDa

40 kDa

31 kDa

22 kDa

Antibody used in WB on Human Liver at 0....

SAMD4A

90 kDa

65 kDa

40 kDa

31 kDa

22 kDa

Rabbit Anti-SAMD4A
Sample Type: Human Fetal Brain
Antibody Concentration: 1 µg/mL

Antibody used in WB on Hum. Fetal Brain ...

SAMD4A

90 kDa

65 kDa

40 kDa

31 kDa

22 kDa

Rabbit Anti-SAMD4A
Sample Type: Human Fetal Liver
Antibody Concentration: 1 µg/mL

Antibody used in WB on Hum. Fetal Liver ...