

Product Datasheet

C15orf27 Antibody (orb1240765)

Description

C15orf27 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Heart at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human C15orf27.

Target

C15orf27

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store C15orf27 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

C15orf27 antibody can be used for detection of C15orf27 by ELISA at 1:312500. C15orf27 antibody can be used for detection of C15orf27 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

39 kDa

Uniprot ID

[Q2M3C6](#)

NCBI

[AAI04954](#)

Dilution Range

C15orf27 antibody can be used for detection of C15orf27 by ELISA at 1:312500. C15orf27 antibody can be used for detection of C15orf27 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.