

Product Datasheet

CBLN4 Antibody (orb1240780)


Description

CBLN4 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CBLN4.
Target	CBLN4
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CBLN4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	CBLN4 antibody can be used for detection of CBLN4 by ELISA at 1:312500. CBLN4 antibody can be used for detection of CBLN4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	Q9NTU7
NCBI	NP_542184
Dilution Range	CBLN4 antibody can be used for detection of CBLN4 by ELISA at 1:312500. CBLN4 antibody can be used for detection of CBLN4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HeLa at 1 u...



Antibody used in IHC on Human cortex at ...



Antibody used in IHC on Human cortex.