



Product Datasheet

HMOX1 antibody (orb12408)

Description	Rabbit polyclonal antibody to HMOX1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human Heme Oxygenase 1 (HO-1/HMOX1) (NP_002124.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	33KDa
Uniprot ID	P09601
Entrez	3162
Dilution Range	WB: 1:200-1:500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western blot analysis of extracts of var...