

Product Datasheet

ACTR1B Antibody (orb1240809)

Description

ACTR1B Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ACTR1B.

Target

ACTR1B

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ACTR1B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ACTR1B antibody can be used for detection of ACTR1B by ELISA at 1:62500. ACTR1B antibody can be used for detection of ACTR1B by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

42 kDa

Uniprot ID

[P42025](#)

NCBI

[NP_005726](#)

Dilution Range


ACTR1B antibody can be used for detection of ACTR1B by ELISA at 1:62500. ACTR1B antibody can be used for detection of ACTR1B by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in WB on Human HepG2 at 0....



Antibody used in WB on Human 721_B at 1 ...



Antibody used in WB on Hum. Fetal Brain ...