

Product Datasheet

CYP1B1 Antibody (orb1240836)

Description

CYP1B1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CYP1B1.

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CYP1B1.

Target

CYP1B1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CYP1B1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CYP1B1 antibody can be used for detection of CYP1B1 by ELISA at 1:62500. CYP1B1 antibody can be used for detection of CYP1B1 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

61 kDa

Uniprot ID

[Q16678](#)

NCBI

[NP_000095](#)

Dilution Range

CYP1B1 antibody can be used for detection of CYP1B1 by ELISA at 1:62500. CYP1B1 antibody can be used for detection of CYP1B1 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90kDa


40kDa

30kDa

20kDa

Antibody used in WB on Human Liver at 0....

CYP1B1



See Immunoblot Data and Customer Feedback for more information

Antibody used in WB on Human lung micros...