

Product Datasheet

ASPH Antibody (orb1240878)

Description

ASPH Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ASPH.

Target

ASPH

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ASPH antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ASPH antibody can be used for detection of ASPH by ELISA at 1:62500. ASPH antibody can be used for detection of ASPH by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

25 kDa

Uniprot ID

[Q12797](#)

NCBI

[NP_064549](#)

Dilution Range

ASPH antibody can be used for detection of ASPH by ELISA at 1:62500. ASPH antibody can be used for detection of ASPH by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90 kDa
65 kDa
40 kDa
31 kDa
22 kDa

Antibody
used in WB
on Human
HepG2 at 0....

ASPH
90
65
40
31
22
Rabbit Anti-ASPH
Sample Type: Human T21_B
Antibody Concentration: 1 µg/mL

Antibody
used in WB
on Human
721_B at 1 ...