

Product Datasheet

BRK1 Antibody (orb1240888)

Description

BRK1 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human C3orf10.

Target

BRK1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store C3orf10 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

C3orf10 antibody can be used for detection of C3orf10 by ELISA at 1:62500. C3orf10 antibody can be used for detection of C3orf10 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

9 kDa

Uniprot ID

[Q8WUW1](#)

NCBI



[NP_060932](#)

Dilution Range

C3orf10 antibody can be used for detection of C3orf10 by ELISA at 1:62500. C3orf10 antibody can be used for detection of C3orf10 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Liver.



Antibody used in WB on Human HepG2 at 2....