


Product Datasheet

RSRC2 Antibody (orb1240891)


Description

RSRC2 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RSRC2.
Target	RSRC2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RSRC2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	RSRC2 antibody can be used for detection of RSRC2 by ELISA at 1:12500. RSRC2 antibody can be used for detection of RSRC2 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	50 kDa, 44 kDa, 22 kDa
Uniprot ID	Q7L4I2
NCBI	NP_075388
Dilution Range	RSRC2 antibody can be used for detection of RSRC2 by ELISA at 1:12500. RSRC2 antibody can be used for detection of RSRC2 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Stomach.



Antibody used in WB on Human Jurkat cell...