


Product Datasheet

MED4 Antibody (orb1240955)


Description

MED4 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MED4.
Target	MED4
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store MED4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	MED4 antibody can be used for detection of MED4 by ELISA at 1:62500. MED4 antibody can be used for detection of MED4 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	30 kDa
Uniprot ID	Q9NPJ6
NCBI	NP_054885
Dilution Range	MED4 antibody can be used for detection of MED4 by ELISA at 1:62500. MED4 antibody can be used for detection of MED4 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HepG2 at 0....



Antibody used in IHC on Human Testis.



Antibody used in IHC on Human Testis.